

LifePack 1000 AUTOMATED EXTERNAL DEFIBRILLATOR (3 LEAD) (CFD PARAMEDICS ONLY)

INDICATIONS

- AED 1000 is to be brought to patients' side for all EMS responses by CFD ALS ambulances.

CONTRAINDICATIONS

- None

EQUIPMENT

- LifePack AED 1000 with case
- 3 Lead cable
- 2 sets of therapy pads
- Electrodes (adult and pediatric)

PROCEDURE

1. Initiation of Patient Care policy B-2.
2. The "3 Lead cable", with electrodes will be applied when indicated i.e. based on patient's history and/or chief complaint.
3. The "therapy pads" will be applied to the patient when indicated for defibrillation.
4. At anytime the LifePack 1000 (3-lead AED) is used, either for monitoring or for treatment, the patient will be immediately transferred to the cardiac monitor/defibrillator, once the patient has been transferred into the ambulance.
5. For known cardiac or respiratory arrest calls, the monitor/defibrillator will be carried into patient's side.

Copyright 2016 Chicago EMS Medical Directors Consortium

Written: 6/11

Reviewed: 6/11

Revised:

MDC Approval: 6/16/11

IDPH Approval: 9/29/11

Implementation: 4/1/12